Shoulder spacer with gentamicin

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Shoulder spacer with gentamicin.

DESCRIPTION:

The shoulder spacer is a temporary implant whose main function is to maintain the articular space during Two Stage Revision Arthroplasties, and to help treating the infection in the prosthesis implant site. The technique consists of two steps: during the first one, the infected prosthesis is removed and the spacer is implanted. During the second step, once the infection has been controlled, the spacer is removed and the definitive prosthesis is implanted.

COMPOSITION:

In order to increase its resistance, the shoulder spacer consists of a central core in AISI 316L stainless Steel, covered with PMMA surgical cement with high Gentamicin concentration.

PROPERTIES OF GENTAMICIN:

Gentamicin is a broad spectrum aminoglycoside antibiotic. It is effective on gram (-) aerobic bacteria, including enterobacteriaceae, Pseudomonas and Haemophilus. It is also effective on Gram positive germs such as staphylococci (Staphylococcus aureus and Staphylococcus epidermidis) including penicillinase producing strains, and has a very limited activity on streptococci. It has no activity on anaerobic bacteria.

MODELS AND GENTAMICIN CONTENT

| MODEL | 40 | 48 |
|---------------------|------|------|
| GENTAMICIN BASE (G) | 1.00 | 1.50 |

WARNING REGARDING

GENTAMICIN:

Patients treated with aminoglycosides must be carefully controlled due to the possible toxicity on kidneys, hearing, and peripheral nervous system which could be associated to the use of this group of antibiotics. The concomitant use with other neurotoxic and/or nephrotoxic medicines such as cisplatin, cephaloridine, kanamycin, amikacin, streptomycin, neomycin, polymyxin B, colistin, paromomycin, vancomycin, tobramycin, viomycin, as well as diuretics like etacrynic acid or furosemide which could enhance its ototoxicity. Old age and dehydration can increase the risk of toxicity in the patient.

PREGNANCY:

Gentamicin falls into the D Category of FDA for medicines: There is positive evidence of human fetal risk based on adverse reaction data from investigational or marketing experience or studies in humans. As pregnant women may benefit from the use of this drug, it's prescription should be evaluated by their treating physicians.

BREASTFEEDING:

Although Gentamicin is poorly excreted into breastmilk it's use should be avoided during this period.

CONTRAINDICATIONS:

The use of PMMA cements with Gentamicin is contraindicated in

patients with a medical history of allergy to any of its components.

USE IN GERIATRICS:

It is known that this drug is excreted in significant amounts by the kidneys and that there is an increase of side effects in patients suffering from kidney insufficiency.

As elderly persons have more possibilities of presenting some degree of kidney insufficiency that might not be manifest in routine laboratory tests, a Creatinine Clearance test can yield more specific information. It is very important to closely monitor the kidney function throughout the treatment.

USE IN PAEDIATRICS:

The use of PMMA cements with or without antibiotics in paediatric patients must be restricted to those exceptional cases in which no other alternative treatments are available.

PRESENTATION:

The external packaging encloses the instructions for use and two successive blisters sealed by a Tyvek® film which contain the shoulder spacer. Sterilization Method: ETO.

INDICATIONS:

The Spacer permits to achieve a high local antibiotic concentration with a low systemic concentration.

Its use is indicated for a short period of time in patients suffering from an infection by Gentamicin sensitive germs. It is recommended to use an appropriate systemic antibiotic therapy together with the spacer. The Spacer is placed as a

hemiarthroplasty, once the original implant has been removed. It is inserted into the humeral canal without cement to allow the release of the antibiotic from the spacer into the surrounding tissues.

In order to properly rehabilitate the patient, it is necessary to evaluate the range of movements within which the stability of the joint is maintained. As the stability of the shoulder joint is closely related to the function of the surrounding muscles, the spacer includes a three hole fin which allows the temporary reinsertion of the rotating cuff muscles providing more stability to the shoulder, a better movement range, less pain and a better quality of life until the definitive prosthesis is re-implanted. The rehabilitation with this definitive prosthesis will also be facilitated because the spacer will help to preserve the rotating cuff and the muscular trophism.

The Spacer is to be kept "in situ" until the surgeon deems it necessary: usually until the infection has been controlled, and it is removed and replaced by a definitive prosthesis. The spacer must not remain implanted beyond 6 months.

CONTRAINDICATIONS:

The use of the Shoulder Spacer is contraindicated in the following cases:

ABSOLUTE:

1. In patients sensitive to Gentamicin or PMMA.
2. In patients suffering from impairments or mental and/or neurologic conditions which prevent them from complying with the post-operative medical indications.
3. In patients who have not

undergone a total shoulder replacement and whose infection is secondary to trauma, septic arthritis, or another surgical procedure.

4. In the presence of Gentamicin resistant pathogens .
5. In the absence of an adequate bone structure allowing the appropriate support of the prosthesis into the proximal humerus or in the glenoid region.
6. In patients with a history of renal pathology, certain diseases in the hearing apparatus and in the neuromuscular end plate.

RELATIVE:

1. In immuno-depressed patients, those suffering from nutrition deficiencies or those whose health is otherwise compromised.
2. In patients with deficiencies in their nervous, vascular or muscular systems.
3. When the etiological agent of the infection could not be identified.
4. In patients with a diagnosis of a systemic secondary or remote infection, or if such infection is liable to develop.
5. In patients with neuromuscular disorders which prevent them from controlling their shoulder joint.
6. When the bone structure is weak, as in an osteoporosis, because it can lead to a migration of the prosthesis, or entail the fracture of the bone.
7. In patients without sufficient bone tissue to permit the insertion and fixation of the prosthesis.

WARNING:

The procedure must be carried out by medical professionals with an appropriate expertise, in an

institutional environment equipped with all the elements required to ensure a normal development of the surgery.

While the shoulder spacer remains implanted it is contraindicated to load weights, carry objects, etc. The surgeon will determine the type and range of movements the patients can make, taking into account the post-operative and clinical condition of the patient.

CAUTION:

- The medical professional in charge must warn the patient that the Spacer may break as a result of certain non-authorized movements or activities.
- As a temporary joint prosthesis the Shoulder Spacer has mechanical limitations which should be properly explained to the patient.
- The decision of implanting a spacer in a pregnant or breastfeeding woman is at the discretion of the medical professional in charge, and requires a risk-benefit evaluation.
- Do not use this product after the expiry date shown in the labelling.
- Single use Product. DO NOT RE-USE. Re-using it may entail mechanical, physical, chemical and/or biological contamination risks.
- After using this product, discard it and dispose of it according to the laws and regulations in force.
- Do not re-sterilize.
- Discard if any of the successive blisters ensuring sterility are found damaged or open.
- Handle with care and store in a cool (below 25°C) and dry place, away

from the light.

- SOLD EXCLUSIVELY TO PROFESSIONALS AND SANITARY INSTITUTIONS.

SURGICAL TECHNIQUE:

During pre-operative planning, it is possible to determine the adequate size of the spacer to be used by comparing radiological images with the X-ray templates supplied to this end.

1. The Shoulder Spacer is implanted using any access way that the surgeon deems appropriate.
2. Remove the original prosthesis. Thoroughly remove all of the cement that was used, as well as any foreign matter found. Resect all tissue that looks devitalized or suspected of being infected. Bear in mind that you should leave as much healthy tissue as possible, in order to permit the re-implantation of the definitive prosthesis once the infection is cured.
3. By means of a metal test spacer, during surgery, verify that the size of spacer selected during pre-operative planning is correct. Finally, implant the spacer that best fits the size of the glene and humerus of the patient.
4. Insert the spacer into the humeral canal and reduce the shoulder as you would in a regular hemiarthroplasty.
5. In some cases it may be necessary to use a small amount of cement on the base of the humeral head in order to ensure a better implant stability.
6. The spherical and polished finish of the humeral head of the spacer prevents the bone stock wear.

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Synergie Ingénierie Médicale S.A.R.L.