## **Knee spacer with gentamicin**

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## Knee spacer with gentamicin

**DESCRIPTION:**

The Knee Spacer is a temporary implant which main function is to preserve the articular space during a Two Stage Revision Arthroplasty. Furthermore it also aids in the treatment of the infection that emerged due to the implant of the knee prosthesis. This procedure consists of a first stage during which the infected prosthesis is removed, and the spacer is implanted. During the second stage, once the patient has finished the specific treatment and no signs of the infection remain, the spacer is removed to be replaced by a definitive new prosthesis.

**KNEE SPACER MODELS:**

RD58 / RD65 / RD72 / RD79 / RI58 / RI65 / RI72 / RI79.

Each measure is offered in two models: Right (D) and Left (I). The surgeon must ensure that the spacer is the size and model right to preserve bone stock.

**COMPOSITION:**

The Knee Spacer is made of bone cement with a high concentration of Gentamicin.

**GENTAMICIN CONTENT:**

| MODEL       | GENTAMICIN BASE (G)(*) |
|-------------|------------------------|
| RD58 / RI58 | 3.2                    |
| RD65 / RI65 | 3.5                    |

|             |     |
|-------------|-----|
| RD72 / RI72 | 5.3 |
| RD79 / RI79 | 5.5 |

(\*) THE INDICATED VALUE MAY VARY IN  $\pm 5\%$

**PRESENTATION:**

The Knee Spacer comes inside a double packaging. The primary packaging consists of two successive blisters sealed with Tyvek lids, to ensure the sterility of the product. The secondary packaging (outer case) contains the primary packaging, the patient's record tab, and the instructions for use. Sterilization Method: ETO

**INDICATIONS:**

The Knee Spacer yields a high local Gentamicin concentration while maintaining low systemic levels. Its use is indicated for a limited time span, in patients suffering from an infection by germs sensitive to gentamicin, requiring a Two Stage Revision Arthroplasty.

A systemic antibiotic therapy should be used in combination while the spacer is implanted.

The Knee Spacer is placed like a regular arthroplasty, after the removal of the original implant.

This knee spacer consists of two independent parts: a tibial plate and a femoral component. The first has a flat base where the femoral component articulates. It is recommended that both components be fixed to the bone by means of a small amount of antibiotic impregnated bone cement. The Knee Spacer maintains the articular space, the length of the affected limb, and the ligament

apparatus of the knee.

The spacer is to be kept in place until it is replaced by the definitive prosthesis when the treating physician deems it appropriate, usually within 6 months of its implantation date.

#### **PROPERTIES OF GENTAMICIN:**

The Gentamicin released by the Knee Spacer is a wide spectrum aminoglycoside antibiotic, with bacteriostatic and bactericidal action against Gram positive and Gram negative germs and anaerobic bacteria.

#### **CONTRAINDICATIONS:**

The following conditions contraindicate the use of the Knee Spacer.

##### **ABSOLUTE:**

- Patients who are allergic or sensitive to Gentamicin or to any of the PMMA components.
- Patients whose mental or neurological condition, make them unable to comply with the postoperative medical indications.
- Patients who have not undergone a total knee replacement and whose infection is secondary to trauma, septic arthritis or any other surgical procedure.
- In the presence of Gentamicin resistant pathogens.
- In the absence of an adequate bone structure which would permit the appropriate support for the prosthesis in the distal femur or in the tibia.

##### **RELATIVE:**

- Immuno-depressed patients, or those with nutritional deficiencies,

or cases where health might be endangered in any other way.

- Patients with deficiencies in their nervous, vascular or muscular systems.
- Cases in which the etiologic agent of the infection could not be determined.
- Patients in which a systemic or remote secondary infection has been diagnosed or when such an infection can be expected.
- Patients suffering from a neuromuscular disorder which prevents them from controlling their knee joint.
- Deficits in the structure of the bone as in patients suffering from osteoporosis, for it might entail a prosthesis migration or a bone fracture.
- Patients without enough bone tissue to permit the insertion and fixation of the prosthesis.
- Patients whose weight or activity levels might result in an early failure in the performance of the spacer.

##### **WARNING:**

- Surgery must be performed by medical professionals with an adequate experience in this procedure, in an appropriate institutional environment equipped with all the necessary elements to ensure the normal development of the procedure.
- Under no circumstances should these devices replace the systemic antibiotic treatment with the appropriate dosage for each case.
- Due to the gentamicin provided

the spacer, the professional in charge must evaluate the risk of using other systemic which might originate ototoxicity or nephrotoxicity.

- As the Knee Spacer is a temporary joint with certain mechanical limitations, it is contraindicated to bear the entire body weight on the affected limb.
- The treating team should prescribe the kind of movements the patient can make, as well as the recommended walking aids such as crutches, canes, etc.
- The patient must be warned that certain movements or activities may entail the breakage or dislocation of the spacer.
- The following auxiliary elements are supplied: X-ray templates and Test Spacers (See Surgical Technique).
- Verify that the test spacers are in appropriate physical and hygienic conditions for their use, in order to prevent a cross-infection or a re-infection.
- Spacers must not be subjected to excessive force during the implantation/insertion (i.e.: striking), as this might produce damages which could reduce the mechanical properties of the device, and not achieve the expected results.
- When removing the spacer special care must be paid not to leave cement residues detached from the spacer.

#### **USE PRECAUTIONS:**

- Single use product.

- Do not re-use. This single use product is designed to be used on a single patient. Re-using it might entail mechanical, physical/chemical, and /or biological contamination risks.
- Handle with care and store cool and dry place.
- Do not use this product after the expiry date shown on the product labeling.
- Discard if any of the successive packaging units which ensure sterility are found damaged or open.
- Do not re-sterilize.
- Due to the gentamicin contents the use of this product is not recommended during pregnancy or likely pregnancy. Its use should only be evaluated if the patient's life is at risk.
- **SOLD EXCLUSIVELY TO MEDICAL PROFESSIONALS AND SANITARY INSTITUTIONS.**
- After removing the spacer, and before the insertion of the definitive prosthesis, a careful cleaning must be performed to eliminate any foreign matter present.

#### **ADVERSE EFFECTS:**

The following adverse effects may be related to the use of spacers:

- Allergic reactions to Gentamicine or to PMMA orthopedic cement.
- Risk of generating ototoxicity or nephrotoxicity due to the gentamicin provided by the spacer, depending on the clinical condition of the patients.
- This product may break if subjected

to excessive weight or non-authorized movements.

- Dislocation or proximal movement due to an inadequate fixation to the bone.

### **SURGICAL TECHNIQUE:**

1. The Knee Spacer can generally be implanted without restrictions though any access way the surgeon deems appropriate.
2. Remove the original prosthesis. Thoroughly remove remaining cement as well as any foreign matter. Resect any tissue that might look devitalized or infected. Bear in mind that it is necessary to leave as much healthy tissue as possible in order to allow the re-implantation of the definitive prosthesis once the infection is cured.
3. At the preoperative planning, by superimposing an X-ray template on the patient's x-ray, it will be possible to select the adequate spacer for the patient. During the insertion procedure, confirm the size selected by means of the reusable Test Spacers to choose the one that best fits the size of the femoral condyles and the tibial plateau of the patient. (Verify that the test spacers are in appropriate physical and hygienic conditions for their use, in order to prevent a cross-infection or a re-infection).
4. The femoral component is to be installed first and must fit adequately. It is recommended to use a small amount of antibiotic impregnated cement can be used to keep it in place. The tibial

component is also implanted using a small amount of antibiotic impregnated cement. A flexion/extension movement should be made before the cement sets so as to permit the self-centering of the device. Remove the excess cement. Reconstruct the extensor apparatus of the knee, maintaining a joint excursion between 0° and 90°.

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Synergie Ingénierie Médicale S.A.R.L.