Hip spacer with gentamicin

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Hip spacer with gentamicin.

DESCRIPTION:

The Hip Spacer is a temporary implant which main function is to preserve the articular space during a Two Stage Revision Arthroplasty. Furthermore it also aids in the treatment of the infection that emerged due to the implant of the hip prosthesis. This procedure consists of a first stage during which the infected prosthesis is removed, and the spacer is implanted. During the second stage, once the patient has finished the specific treatment and no signs of the infection remain, the spacer is removed to be replaced by a definitive new prosthesis.

HIP SPACERS MODELS

MÜLLER TYPE: 48M / 56M.

CHARNLEY TYPE: 40C / 48C / 56C / 40C XL / 48C XL / 56C XL.

COMPOSITION:

The Hip Spacer has an implantable AISI 316 L stainless steel core -which increases its mechanical resistance, covered with PMMA bone cement with a high Gentamicin concentration.

GENTAMICIN CONTENT

MODEL	GENTAMICIN BASE (G)(*)
40C	1.60
48C	2.50
48M	2.70
56C	3.70
56M	3.80
40C XL	2.10
48C XL	3.00

56C XL

4.20

(*) THE INDICATED VALUE MAY VARY IN $\pm 5\%$

PRESENTATION:

The Hip Spacer is contained in a double barrier tray. Trays are sealed with Tyvek® in order to ensure the sterility of the device.

Trays are protected by a cardboard box also containing the Instructions for Use. Sterilization Method: ETO

INDICATIONS:

Its use is indicated during a limited time-span, in patients suffering from an infection by Gentamicin sensitive germs, requiring a Two Stage Revision Arthroplasty.

A systemic antibiotic therapy should also be prescribed while the spacer remains implanted.

The Hip Spacer is implanted after the removal of the infected implant, as a regular Hemiarthroplasty. The spacer is inserted into the femoral canal without using any cement in order to allow a better release of the antibiotic to the surrounding tissues. A handmade ring of PMMA bone cement may be added to the base of the spacer's neck in order to increase stability.

The Hip Spacer maintains the joint space and the length of the affected limb, so preserving the abduction and stabilization apparatus of the hip. The spacer is to be kept in place until it is replaced by the definitive prosthesis when the treating physician deems it appropriate, usually within 6 months of its implantation date.

PROPERTIES OF GENTAMICIN:

The Gentamicin released by the Hip Spacer is a wide spectrum aminoglycoside antibiotic, with bacteriostatic and bacteriocidal action against Gram positive and Gram

negative germs, and anaerobic bacteria.

CONTRAINDICATIONS:

The following conditions contraindicate the use of the Hip Spacer:

ABSOLUTE:

1. Patients who are allergic or sensitive to Gentamicin or to any of the PMMA components.
2. Patients whose mental or neurological condition, make them unable to comply with the postoperative medical indications.
3. Patients whose mental or neurological condition, make them unable to comply with the postoperative medical indications.
4. Patients who have not undergone a total hip replacement and whose infection is secondary to trauma, septic arthritis or any other surgical procedure.
5. Lack of bony structure capable of providing adequate support to the spacer's rod and / or head.

RELATIVE:

1. In the case of immunodepressed patients, those with nutritional deficiencies, or those whose health might be endangered in any other way.
2. In patients with deficiencies in their nervous, vascular or muscular systems.
3. In those cases where the etiologic agent of the infection could not be determined.
4. In patients being diagnosed with a systemic or remote secondary infection, or when such an infection is liable to develop.
5. In patients suffering from a neuromuscular disorder which prevents them from controlling the

hip joint.

6. When the bone is weak, such as in patients suffering from osteoporosis, for there might be a prosthesis migration or a bone fracture.
7. In patients without sufficient bone tissue to permit the insertion and fixation of the prosthesis.
8. In patients whose weight or activity levels might entail an early failure in the performance of the spacer.

WARNING

- Surgery must be performed by medical professionals with an adequate experience in this procedure, in an appropriate institutional environment equipped with all the necessary elements to ensure the normal development of the procedure.
- Under no circumstances should these devices replace the systemic antibiotic treatment with the appropriate dosage for each case.
- Due to the gentamicin provided the spacer, the professional in charge must evaluate the risk of using other systemic which might originate ototoxicity or nephrotoxicity.
- As the Hip Spacer is a temporary articulation, with certain mechanical limitations, it is contraindicated to bear body weight on it.
- The professional team should indicate the kind of movements the patient can make during rehab, as well as prescribing and supervising the use of weight bearing devices such as crutches, canes, etc.
- The patient must be warned that certain movements or activities may entail the breakage or dislocation of the spacer.
- The following auxiliary elements are

supplied: X-ray templates and Test Spacers (See Surgical Technique).

- Verify that the test spacers are in appropriate physical and hygienic conditions for their use, in order to prevent a cross-infection or a re-infection.
- Spacers must not be subjected to excessive force during the implantation/insertion (i.e.: striking), as this might produce damages which could reduce the mechanical properties of the device, and not achieve the expected results.
- When removing the spacer special care must be paid not to leave cement residues detached from the spacer.

PRECAUTIONS

- Single use product.
- Do not re-use. This product is designed for a single use on a single patient. Re-using it might entail mechanical, physical/chemical, and / or biological contamination risks.
- Handle with care and store in a dark and dry place.
- Do not use this product after the expiry date shown on the product labeling.
- Do not re-sterilize.
- Discard if any of the double trays which ensure sterility are damaged or open.
- Due to the gentamicin contents the use of this product is not recommended during pregnancy or likely pregnancy. Its use should only be evaluated if the patient's life is at risk.
- **SOLD EXCLUSIVELY TO MEDICAL PROFESSIONALS AND SANITARY INSTITUTIONS.**
- After removing the spacer, and before the insertion of the definitive prosthesis, a careful cleaning must be performed to eliminate any foreign

matter present.

ADVERSE EFFECTS

The following adverse effects may be related to the use of spacers:

- Allergic reactions to Gentamicin or to PMMA orthopedic cement.
- Risk of generating ototoxicity or nephrotoxicity due to the gentamicin provided by the spacer, depending on the clinical condition of the patients.
- This product may break if subjected to excessive weight or non-authorized movements.
- Dislocation or proximal movement due to an inadequate fixation to the bone.

SURGICAL TECHNIQUE

1. The Hip Spacer can generally be implanted without restrictions, though any access way the surgeon deems appropriate.
2. Remove the original prosthesis. Thoroughly remove all cement leftovers as well as any foreign body. Resect any tissue that might look devitalized or infected, considering as much healthy tissue as possible should be left in place in order to allow the re-implantation of the definitive prosthesis once the infection is cured.
3. During the preoperative planning, the adequate size of spacer can be determined by means of x-ray templates. During the implantation procedure, confirm the size that better suits the patient by means of the Fitting Spacers.
4. Insert the spacer into the femoral canal and reduce the hip as in a traditional hemiarthroplasty. A well preserved acetabulum can provide

- a direct support to the Hip Spacer.
- In case of a bone stock deficiency, the support of the spacer on the remaining acetabulum should be evaluated, while the patient's mobility is accordingly restricted.
 - The perfectly spherical and polished finish of the femoral head of the spacer prevents the bone stock wear.

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Synergie Ingénierie Médicale S.A.R.L.