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Opening wedge high tibial osteotomy plate

Fixation monobloc™ Single unit fixation™

OSTÉOTOMIE
OSTEOTOMY

TECHNIQUE OPERATOIRE
SURGICAL TECHNIQUE

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The SURFIX® osteotomy plates (PTOD and PTOG) are used for fixation of opening wedge high tibial osteotomy and for the treatment of osteoarthritis of the tibio-femoral compartment with genu-varum.

These plates may be also used for:


- internal fixations of fractures of the medial tibial plateau,
- closing internal wedge tibial osteotomy of varisation.

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The cutting guide for opening wedge osteotomy

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The SURFIX® opening cutting guide provides a precise and guided osteotomy line and allows to measure the depth of the cut.

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
The SURFIX® distractor allows opening at the correct angulation, maintaining this correction while filling up the gap and the fixation of the SURFIX® plate for opening wedge osteotomy.

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The surgical technique

SURFIX®, the manufacturer of this device, does not practise medicine and does not recommend any particular surgical technique for a given patient. Only the surgeon fitting this device is responsible for deciding on and using the appropriate techniques for each patient.

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
Preoperative planning is performed with long leg standing frontal plane X-rays. Measurement is done using the anatomic axis method, measuring the angle between the anatomical axis of the femur and the tibia (HKA angle) **a**.

The radiological technique used is to be as accurate and reproducible, according to a protocol described and controlled by the surgeon **b**.

The aim of the surgery is to accomplish angular correction of 3° to 6° of valgus post operatively, so as to increase the chance for consolidation and for durable results.

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The patient is placed supine on the operation table. A sang bag is placed under the involved hip. A leg holder is applied to the proximal thigh, maintaining the knee in a semi flexed position. The knee is thus positioned with the patella being oriented in a strict frontal plane. A tourniquet is applied to the proximal thigh. The C-arm image intensifier is placed in such way that the femorotibial joint is well seen in the frontal plane.

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3


Approach and exposure in the osteotomy area

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Vertical metaphyseal anterior-interior incision, slightly sinuous or arciform with posterior concavity, mid distance between the anterior and posterior borders of the medial tibia cortex, not extending beyond the level of the joint.

Section of the subjacent aponeurotic plane along the same path and detachment in one piece, using a scalpel or a sharp bone chisel, of the anterior-superior part of the

insertions of the pes anserinus and the foot of the superficial medial collateral ligament, from front to back; thus medial tibial cortex of the tibial metaphysis is exposed over a height of 4-5 cm.

Release of the posterior face of the metaphysis with a blunt periosteal elevator (e.g. Lambotte) which will be left in place, or replaced by a Hohmann or straight Merle D'Aubigné retractor.

Freeing of the patella tendon at its insertion, and the subtendinous bony area in the direction of the outer facet of the metaphysis.

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
Using a motorised tool, placement of one or two K-wires (diam. 1.8 mm) **c**, from inside to outside, with the help of the SURFIX® cutting guide for opening wedge osteotomy positioned previously:

- external tip inserted as far as the bone, at approximately 15 mm under the joint **d**, marked by an IM needle inserted in the latter.

- stabilised blade guide in contact with the medial cortex **e**. The most posterior K-wire must be parallel to the plane of the posterior cortex of the metaphysis. The level and obliquity of the K-wires are checked by C-arm.

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
The saw blade must be wide and rigid, to avoid the risk of a change of plane. It is inserted in the slot of the SURFIX® cutting guide for opening wedge osteotomy **f**, and will cut nearly the whole of the tibial metaphysis, taking care at the front not to damage the foot of the patella tendon, and at the back not to cut the posterior cortex. To prevent the risk of posterior escape of the blade the xyster (or retractor) is always in place to protect soft tissues, a lateral cortical hinge of 3-5 mm must be retained. The width of the metaphysis can be read directly on the graduated section of the SURFIX® cutting guide for opening wedge osteotomy **g**.

The depth of the cut can be measured by reading the measurement on a graduated saw blade after deducting the thickness of the blade guide (14 mm).

This allows to:

- assess the thickness of the osteotomy hinge,
- select the ruler for the SURFIX® distractor.

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
The cutting guide for opening wedge osteotomy is dismantled and removed, leaving the pins in place.

Use a very sharp osteotome, 25-30 mm wide, allowing the osteotome to penetrate into the metaphysis little by little, respecting the orientation of the guide **h**.

The osteotomy must conserve a lateral cortical hinge: tactile and auditory perception of cortical contact or visual check by fluoroscopy ensure that it is maintained.

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The depth of the osteotomy site may be measured either by using the cutting guide for opening wedge osteotomy, or by using a K-wire and rule.

The distractive forceps are fitted with one of the two small rulers whose face corresponds best to the depth of the osteotomy site **i**.

The angle of the opening may be read directly from the small ruler **j**. (With reference to "Open wedge tibial osteotomy: combined coronal and sagittal correction.

Ph. Hemigou - The knee 9(2002)15-20.

Gradually move the leg segment into valgus, the lateral side of the thigh immobilised by a support positioned previously (or with the assistance of manual counter-support), to slightly open the osteotomy site.

Insert the tips of the SURFIX® distractor horizontally, on the corner between the medial and posterior cortex **k**, then open gently the osteotomy site with the forceps until the desired correction angle is reached **l**.

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This is at the surgeon's discretion, using either an iliac autograft or synthetic bone substitute.

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7.a - By autograft

The osteotomy opening is measured in millimetres, at the front and at the back. It is usually slightly asymmetrical, wider at the back than at the front (by exaggerating this difference by placing the knee in hyper-extension, a slight pre-operative flexum of arthritic origin may be corrected). Extraction of 2 triangular cortical-cancellous grafts on the homolateral iliac crest of the desired thicknesses and a few cancellous shavings using a curette or gouge.

Insertion of the narrower graft in the anterior side of the osteotomy, removal of the distractor, and insertion of the wider graft at the back, temporarily increasing the opening to facilitate insertion. The cancellous bone shavings are packed between the grafts to fill any residual gaps.

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7.b - By synthetic bone substitute


SURFIX® offers a range of ATAO® anatomic wedges composed of 60% hydroxyapatite and 40% tricalcium phosphate with totally interconnected porosity.

The ATAO® substitute must be chosen using the SURFIX® trial for bone substitute. This trial for bone substitute is presented into the osteotomy site and the substitute size to be selected may be read on it **m**. The ATAO® substitute is then inserted using the special SURFIX® bone substitute holder **n**.

The bone substitute is put in place at the surgeon's discretion:

- either before placement of the plate, keeping the distractor in position to hold the opening. The distractor allows the opening to be increased slightly to facilitate insertion of the substitute,
- or after placing the plate (limited access due to the plate).

SURFIX® also offers ATAO® granulate of the same composition as the substitutes, which provide rehabilitating throughout the osteotomy site.

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
Position a PTO SURFIX® plate (D for right or G for left) previously fitted with two 3.5 mm diameter drill guides screwed **o** :

- on the proximal antero superior hole **p**,
- on the inferior distal threaded hole. The plate is usually applied easily on the tibial cortex **q**.

If necessary, it can be remodelled using SURFIX® plate benders.

Screw placement requires certain precautions to ensure that the screws and plate form a perfect single unit, which is the very principle of the SURFIX® process:

- respect the coaxiality between each screw and the corresponding hole in the plate (hence the need to use the drill guides),
- remove any soft tissue or bony fragment between the plate and the screws and at the site of the lock screws **r**,
- Insert the lock-screw just after the screw **s** and lock it tightly **t**,
- fit the screws in a precise order **u**:
1 - cancellous anterior epiphyseal screw,
2 - inferior cortical screw,
3 - intermediate cortical screw,
4 - lastly, cancellous posterior epiphyseal screw.

The 6.5 mm diameter cancellous bone screws may be used on the metaphyseal-diaphyseal part. A 4.5 mm diameter drill guide should then be screwed on the diaphyseal holes and the drill 4.5 mm diameter should be used. It could be useful to use the reamer of the screwdriver for preparing the space for the head of the screw of the cortex bone.

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Collapse of the lateral cortical hinge :

More often than not this does not implicate any additional difficulties, since the high stability provided by the SURFIX® process easily makes up for the potential instability caused by a loss of cortical continuity.

In the case of a large osteotomy opening, if the breakage is total, lateral displacement may be observed at the osteotomy site. This should be corrected by a reverse reduction manoeuvre and stabilised using a single ascending SURFIX® screw transfixing the osteotomy site, avoiding the risks of articular breakage. The SURFIX® plate osteosynthesis is then performed as before.

Fracture of the lateral tibial plateau :

If the osteotomy has not been brought sufficiently outwards, the valgus correction manoeuvre may lead to a vertical split line ending in the lateral tibial plateau, delimiting a fragment of variable size.

The important thing is to notice this in good time. The osteotomy should be closed and the cutting plane resumed with the osteotome on the lateral side, then the osteotomy should be re-opened carefully leaving the osteotome in place, which minimises the risks of displacement of the fracture line.

In the case of displacement or an intra-articular insert deemed not to be tolerable, a percutaneous osteosynthesis of the lateral tibial plateau may be performed using a 6.5 mm diameter SURFIX® cancellous bone screw placed parallel to the plane of the plateau, by lateral approach.

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After haemostasis when the tourniquet is removed.


Re-insertion of the fibre mass corresponding to the termination of the pes anserinus, above the plate, using separate Vicryl 0 sutures, with suction drain. Surface closure in 2 planes according to the surgeon's customary practice.

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
The high stability of the assembly:


- avoids post-operative immobilisation,
- facilitates passive and active mobilisation from the first few days up to 90° flexion without difficulty,
- allows resumption of walking with partial weight bearing assisted by crutches or a walking frame, as soon as the drain is removed.
- The maintenance of a removable posterior splint can be necessary at the heavy patients and little mobiles.

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- Drill guide diam 3.5 and 4.5 mm **1**
- Length gauge **2**
- Drill diam 3.5 mm **3**
- Drill diam 4.5 mm **4**
- Plates benders **5**
- Hexagonal screwdriver diam 3.5 mm **6**

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- Cutting guide for opening wedge osteotomy **7**
- K-wire diam. 1.8 mm **8**
- Trial for bone substitute **9**
- Distractor **10**
- Bone substitute holder **11**
- Trial plate **12**

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Instructions for use SURFIX® Osteosynthesis System

In accordance with ECC directive 93/42 related to medical devices, this product must be handled and/or implanted by well-trained, qualified persons, aware of these directions for use.

1. Packaging and traceability

- Most SURFIX® implants are supplied sterile. The sterilization method is specified on the packaging. Components sterilized by radiation are exposed to a minimum of 25 kGy of g irradiation.
- References and batch numbers are engraved on each implant to facilitate traceability.
- For sterile implants, several copies of labels are provided for the patient's record in each package.

2. Warnings and precautions for handling

- SURFIX® implants must be stored in a dry, cool place away from dust and UV radiation.
- Check that the packaging is intact before using sterile implants. Implants from damaged packages should not be implanted and are considered to be non-sterile. It is recommended that the implant be re-run through the sterilization procedure (see §3).
- The sterility expiration date is limited to five years from the sterilization date (see label).
- Check that the implant matches the specifications mentioned on the package (reference and size) when you unpack it.
- Do not put the implant in contact with objects or substances that might alter its surface or chemical composition. Perform a visual check of each implant before use to detect damage (such as warping, scratches, etc...). In such cases, do not use the implant.
- Inner packaging should be handled under sterile conditions.

3. Recommendation for (re)sterilization

SURFIX® brand products are sold sterile. Products sold sterile, whose sterility is no longer guaranteed (see §2), are considered to be non-sterile and must be (re)sterilized.

- Recommendations for decontamination and cleaning

SURFIX® brand products considered as non-sterile must be unpacked and decontaminated before implanting, in compliance with current regulations. Any stained and/or damaged instruments should be removed from the operating area. Sterilization is not a substitute for cleaning and decontamination.

- Recommendations for (re)sterilization

Re-sterilization is only allowed for non implanted products. Remove delivery packaging in compliance with current regulations to (re)sterilize non-sterile products. SURFIX's osteosynthesis implants are recommended to be sterilized by the steam autoclaving procedure regularly used in the hospital. The suggested method is based on HIMA or AORN recommended practices. The implants can be sterilized several times in the same conditions:

270°F/132°C - 45 minutes	or	298.8°F/134°C - 18 minutes
--------------------------	----	----------------------------

Warning: SURFIX cannot guarantee sterility for implants that have been cleaned or (re)sterilized by the purchaser or user.

4. Product description

SURFIX® osteosynthesis plates are made of 316 L stainless steel for implants in compliance with standard NF ISO 5832-1 and ASTM F138 & F139. SURFIX® osteosynthesis screws and lock screws are made of a TA6V titanium alloy for implants in compliance with standard NF ISO 5832-3 and ASTM F136.

5. SURFIX® fixation system

The SURFIX® system creates extemporaneously a single implant/screw unit fixed into the bone. The osteosynthesis screws must be driven into the bone through the holes in the plate. The system is locked by means of lock screws drilled into the threaded lipped socket at the top of each hole, thus blocking each screw head. The system loses some of its mechanical qualities if not properly locked and much more if it is not locked at all. The SURFIX® plate will then have the same mechanical features as those of a standard screwed plate. The synthesis of the plate onto the bone with the SURFIX® system is not only secured by a tight union between the plate and the bone - as is the case with an ordinary screwed plate - but also because of the good adjustment of the screws into the bone. The screws must be adjusted into the bone as accurately as possible. The screws must be driven through both cortical areas in the case of a synthesis in the metaphyso-diaphysary area.

6. Combined implants

A SURFIX® osteosynthesis implant consists of a plate, as many screws as threaded lipped sockets on the plate and as many lock screws as implanted screws. All these components must bear the SURFIX® trademark. All SURFIX® implants must be fitted with instruments and accessories from the same manufacturer and designed for that purpose. The choice of the best-fitting SURFIX® screw-plate combination is the responsibility of the surgeon.

7. Device performance

- Good stability and seal-like fixation on the two fragments to be synthesised.
- Consolidation of the focus without any compression of the plate on the bone.
- Suppression of the lever effect where the screw and the plate are in contact.

8. Precautions to be taken for fixation

- The SURFIX® system is fully effective only if the screw is perfectly tightened to the implant. The following steps are essential:
- Maintain the performance of the locking system by using only suitable instruments. The bottom of the lipped socket and its thread must be protected so as to avoid damaging the fixation. The oval-shaping of the lipped sockets, when bending the plates into shape, is particularly risky. Some plates have thinner and smaller bendable sections so that screws can be better adjusted. The bending of these sections is carried out with SURFIX® pliers designed for this purpose. It is imperative that the bending be implemented between two consecutive lipped sockets. If this is not the case, the intermediary socket might become oval-shaped. Bend the plates only once and not excessively. The screws must not touch each other once the lipped sockets have been shaped. Bend the plate slightly between two sockets outside the bendable sections to adapt the shape of a plate onto the bone. Plate benders should be used, provided the sockets are not twisted, as mentioned above.
- Drive the screw co-axially through the threaded lipped socket, by using SURFIX® drill guides, which can be fitted into any SURFIX® plates. The drill guide and the drill must have the same diameter.
- The screw must remain co-axial until it has been locked by the lock screw. Consequently, it is necessary to:
- Make sure that the screw can be easily driven into the bone until its collar has reached the bottom of its hole. It may be necessary to widen the top of the hole in the bone. Be careful not to damage the threaded lipped socket, if this is the case.
- Stop screwing as soon as the screw head has reached the bottom of its hole. Further screwing could draw the plate against the bone by exerting a pull-in effect if the implant is not in contact with the bone. This could change the axis of both the screw and the socket.
- Implement each fixation in succession, from the guided drilling to the tightening of the lock screw. The tightening of the screw by the lock screw should always be carried out immediately after the positioning of the screw and before implementing any other operation that could alter the position of the plate against the bone and, consequently, against the screw.
- Remove any foreign bodies that may be found between the implant, the screw and the lock screw.
- Tighten the lock screw.

9. Safety instructions for an osteosynthesis

A strain is created on the screws if the system is fitted to a non-reduced focus and is actually used to reduce it. This can deteriorate the bone and even result in a spiroid fracture or weaken and even break the screws. Therefore, it is strongly recommended to reduce the fracture first, and then fix it with a SURFIX® implant, as described above.

10. Indications

The indications for surgery of the whole range of SURFIX® implants are described in the attached document Indications of these directions for use and in the SURFIX leaflet®.

11. Contraindications

- Acute or chronic, local or systemic infections.
- Lack of musculo-cutaneous cover, severe vascular deficiency touching the focus.
- Bone deterioration preventing a good fixation of the screws in the bone.
- Muscular deficit, neurological deficiency or behavioural disorders, which could submit the osteosynthesis to abnormal mechanical strains.

12. Factors likely to compromise the success of implantation

- Serious vascular deterioration, bone devitalization, severe osteoporosis, loss of bone substance.
- Deformity or severe traumatism with loss of bone substance or soft parts.
- Local bone tumor.
- Systemic, metabolic or genetic disorders.
- Infectious diseases.
- Drug and/or alcohol and/or smoke and/or medications addiction and/or abuse.
- Obesity.
- Intense physical activity (e.g., competitive sports or strenuous work).

The stability of SURFIX® osteosynthesis lessens the painful effect and may therefore induce the patient to move carelessly. The surgeon must warn the patient of the risks of such careless behavior. The risks in connection with these injuries include, for example, torsional stress of the plate and breaking of the screw. These risks would be caused by the abnormally high strain sustained by the osteosynthesis material, particularly with regard to synthesis of highly comminuted fractures with loss of bony substance.

It is important to warn patients having such contraindications about the influence these factors may have on the successful outcome of the operation. Patients should be given advice about the measures to be taken to reduce the effects of these contraindications.

13. Adverse reactions

- The most typical and common adverse reactions following the fitting of osteosynthesis implants are the following:
- Delayed consolidation, pseudoarthrosis.
- Rupture or deformity of all or part of the implant.
- Infection, hematoma, venous thrombosis, pulmonary embolism, cardiovascular problems.


14. Liability

SURFIX TECHNOLOGIES shall not be liable for any accidental or consequential loss, damage or expense, directly or indirectly arising from the use of this product. SURFIX TECHNOLOGIES neither assumes nor authorizes any other person to assume for it any other or additional liability or responsibility in connection with the product. SURFIX TECHNOLOGIES intends that these devices should be used only by physicians having received appropriate training in orthopedic and traumatologic surgery techniques.


IMPORTANT: Please contact your supplier for further information concerning the use of this product.

WARNING: Federal law (USA) restricts this device to sale by or on the order of a physician.


WARNING: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

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

Chamfer the head of the screw with the screwdriver.
Ensure that the threaded hole is not damaged when performing the chamfering.

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
- Insert the screw into the prepared cavity until it reaches the end.
Clean the threaded hole before and after introducing the screw.
Maintain co-axiality between the screw and the threaded hole.

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Insert the lock-screw in the peg which is designed for this purpose.
The screw and lock-screw must be inserted in the same phase of the procedure.

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Lock the lock-screw tightly in its cavity.

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
- All the products are manufactured and reference within the frame of the standards in force.
- Implantation procedure are described in the surgical technique.
- Non-contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- WARNING : Federal law (USA) restricts this device to sell by or on the order of a physician.

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