## EN COLINK® Standard Continuous Loop Instructions for Use

**DEVICE DESCRIPTION:** Teknimed COLINK® Standard Continuous Loop is comprised of a braided ultra-high molecular weight polyethylene (UHMWPE) continuous loop combine with a titanium (Ti-6Al-4V ELI per ASTM F136) plate. Additional non-absorbable sutures consisting of UHMWPE, nylon, polyester, or polypropylene are looped through the titanium plate to aide in passing the plate through the intended void. The UHMWPE is available undyed (white) or with trace filaments of black, blue, or green color suture. The sutures are available uncanted and coated with uniform polybutylene adipate or silicone.

**INDICATIONS:** The Teknimed COLINK® Standard Continuous Loop is intended for use in the fixation of bone and soft tissue in orthopedic procedures requiring ligament or tendon reconstruction.

**CONTRAINDICATIONS:**

- Do not use for surgeries other than those indicated
- Hypersensitivity to the implant material. In patients with suspected hypersensitivity, appropriate tests should be conducted to rule out sensitivity to the implant materials prior to implantation.
- Inadequate bone quality
- Inadequate bone quantity at implantation site.
- Active infection, or previous infections which may retard healing.
- Blood supply limitation
- Patients unable or unwilling to follow directions during the healing process.

**ADVERSE EFFECTS:** Adverse effects associated with the use of this device include: wound dehiscence, callus formation when prolonged contact with salt solutions occur, enhanced bacterial infectivity, minimal acute inflammatory reaction, transitory local irritation, and other reactions to device materials.

**WARNINGS:**

- The device is provided sterile for single use only. Reuse of device may lead to device failure, or patient exposure to harmful contamination leading to infection, serious injury, or death.
- Do not Resterilize.
- Discard if device packaging is opened or damaged, discard unused suture.
- This device should be used by or on the order of a physician. User should be familiar with surgical procedure, techniques, proper implant sizing, and instructions for use prior to use. Appropriate instrumentation should be used for device implantation.
- Preoperatively the instructions for use and device limitations should be communicated to the patient.
- The fixation of the device should be considered temporary postoperatively during the healing process. Postoperative regimen prescribed by the physician should be followed; failure could compromise the device and/or outcomes of the procedure.
- The Teknimed COLINK® Standard Continuous Loop device has not been evaluated for MR related heating or migration. Use of this device in an MR environment may cause device migration, MR related heating, and the MR image may be compromised close to the position of the device.
- A decision to remove the device should take into consideration risks associated with additional procedures, and postoperative management should follow.

**PRECAUTIONS:**

- Care should be taken to ensure proper device placement and to avoid damaging of the bone or device.
- Titanium plate must be properly oriented with the drill channel when passed through bone for correct placement and to avoid damage of bone.
- In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Do not adjust channel size with the device implanted.

**STORAGE CONDITIONS:** The device should not be used past the expiration date and should be stored at ambient conditions.

**HOW SUPPLIED:** Suture supplied meet United States Pharmacopeia (USP) requirements for non-absorbable suture except for diameter, suture dye may include D&C Blue No. 6, D&C Green No. 6, and Logwood Black. The device is sterilized by EO gas, and is provided sterile for single use.

**INSTRUCTIONS FOR USE:**

- Prepare and measure the femoral channel. Use 4.5mm - 7.0mm drill.
- Measure the length of the graft.
- Determine the intra-articular distance from the graft location of the tibia to the femoral socket. Subtract the graft length from the intra-articular distance, and then divide by two to determine the graft insertion length. This ensures the graft insertion is equal on the tibial and femoral sides of the construct.
- The COLINK® Standard size is determined by subtracting the insertion length from the total femoral channel length.
- Prepare the graft loop construct.
- Use the additional sutures to pass the titanium plate parallel through the femoral channel; care should be taken to prevent damage of the femoral channel. Do not advance the titanium plate past the femoral cortex into the overlying tissue.
- Use second suture to flip and rotate titanium plate.
- Tension the tibial side of the construct.
- Confirm final placement on the femoral cortex using X-ray or fluoroscopy.

US Federal Law Restricts this device to sale by or on the order of a medical physician.

		
DO NOT RESTERILIZE	SINGLE USE ONLY	ATTENTION: READ INSTRUCTIONS
	ETHYLENE OXIDE STERILIZATION	
	Manufactured By: Riverpoint Medical 825 NE 25th Ave Portland, OR 97232 Ph. (503) 517-8001 Fax. (503) 517-8002	
	Emergo Europe Westervoortseedijk 60 6827 AT Arnhem The Netherlands	

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