Synicem Clous iSet

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
Instruments for femoral or tibial nails

PART # / REF / (QTY)

DESCRIPTION / USE

01 REF 113166 (1)

GUIDE BODY WITH TWO ADJUSTMENT KNOBS.

(precision guide for the insertion of the FEMORAL/TIBIAL implant and its appropriate proximal locking)

This guide has marks on the front and back faces with the following inscription: (FEMORAL RIGHT / TIBIAL LEFT) or (FEMORAL LEFT / TIBIAL RIGHT).

These marks must be placed so that they are visible at all times during surgery, and correspond with the bone being treated. On the lateral face of the guide there are two holes marked T (tibial) and F (femoral), which permit the manual introduction of a drill bit (supplied in the box) in order to verify the position of the proximal end of the nail being used.

Important: As the diameters of the coupling and the coupling screws are different for femoral and tibial nails, it will not be possible to assemble the device if the guide is not properly set up: Check that the configuration of the guide (coupling and coupling screw) agrees with the bone being treated, Femur or Tibia.

02 REF 113107 (1)

FEMORAL COUPLING SCREW.

03 REF 113168 (1)

FEMORAL COUPLING.

04 REF 113170 (2)

COUPLING NUT (FEMORAL/TIBIAL).

05 REF 113167 (1)

TIBIAL COUPLING SCREW.

06 REF 113169 (1)

TIBIAL COUPLING.

08 REF 113138 (2)

DRILL BITS Ø4MM X 280MM.

For locking through the drill bit guide

09 REF 010530 (2)

DRILL BITS Ø4MM X 170MM.

For distal manual locking.

11 REF 113137 (1)

DEPTH GAUGE FOR LOCKING SCREWS.

This item must be inserted through the locking screw sleeve #23, using the end of this sleeve to read the depth.

12 REF 113128 (1)

DOUBLE TIP SCREWDRIVER.

It can be used both for the locking screws and the nail plug.

13 REF 113120 (1)


WEIGHT FOR SLIDING HAMMER.

14 REF 113116 (1)

SLIDING HAMMER ROD.

15 REF 113133 (1)


SLIDING HAMMER HANDLE.

Note: Items 13, 14 Y 15 are used to transmit forces to the axis of the nail so as to facilitate its insertion.

Warning: The instruments must never be struck directly under any circumstances.


16 REF 113266 (1)

TRIANGULAR TIP AWL.

It is used to initiate the access into the medullary canal.


19 REF 113108 (1)

GUIDE HANDLE.

(Optional) It is used to facilitate the insertion of the implant.

20 REF 113271 (1)

COMBINED #19 WRENCH.

To tighten and loosen: a) the connecting screw, b) the sliding hammer rod, c) the cube.

21 REF 113127 (1)

PUNCH FOR SOFT TISSUES.

For the insertion of the locking screws sleeves.

22 REF 113119 (1)

GUIDES FOR DRILL BITS.

They are screwed into the locking sleeves.

23 REF 113117 (1)

LOCKING SCREWS SLEEVES.

These sleeves are placed and fastened into the body of the guide.

27 REF 010481 (1)

STEINNMAN Ø 3.5MM.

Optional for distal manual locking.

28 REF 113263 (1)

RADIOPAQUE RULER FOR FEMORAL NAILS.

For measurements under X-rays.

29 REF 113264 (1)

RADIOPAQUE RULER FOR TIBIAL NAILS.

For measurements under X-rays.

Accessory box for tibial nail (optional, on demand, not included with the Syncem Cloud iSet)

For Ø8mm tibial nails, to be used with the standard *Syncem Cloud iSet* box.
1 Steinmann, Ø 3.0mm / 2 drill bit guides of Ø3.0mm / 1 Screw box with locking screws of Ø3.9mm and various lengths / (1) Punch for soft tissues / (1) Locking screws sleeves.

Indications:

These *Syncem Cloud iSet* instruments have been exclusively designed to be used for the implantation of Syncem endomedullary spacer Nails for long bones (femur/tibia), following the indication for use of these implants.

Presentation:

The instrument set for the implantation of *Syncem Cloud iSet* femoral or tibial Endomedullary nails is supplied in perforated trays with a specific support and identification for each part.
The trays are packed inside an instrument box which is also multi-perforated and is provided with safety locks. The corresponding assembling, disassembling and post-surgery reprocessing and cleaning instructions are also included.
This set must be sterilized by autoclave (humid heat)

Description:

The instruments for the implantation of a SYNIMED femoral or tibial nail have been designed to be used with a specific range of products under this brand name. Both the box and the tray have been designed for the packaging and protection of the instruments and these must therefore be kept in their original places in the tray. Once the surgical procedure is over, the instruments must be decontaminated and cleaned before putting them back in their corresponding supports. (see illustration/picture).
These products have been manufactured with standardized materials and/or those with an acknowledged biocompatibility.
The user must not modify these elements in any way, either partially or totally.
The maintenance, repair or modification of these products must be performed exclusively by Synergie Ingénierie Médicale SARL or by a firm duly authorized by SYNIMED.

Warning:

None of the devices mentioned herewith are implantable.

All the items contained in this box have been engraved to facilitate their immediate visual identification.

Important:

- The instruments for *Synimed* femoral and tibial nails are

to be used in a surgical environment with all the aseptic conditions described in these instructions for use, in compliance with ISO 8828, and both the surgeons and the surgery auxiliary staff must be familiar with the surgical techniques to be used.

- Both the instruments and the auxiliary items have been exclusively designed to be used in the implantation and removal of *Synimed* tibial and femoral nails. They must not be used with products from other manufacturers under any circumstances. Ignoring this warning might entail the breakage or premature deterioration of the instruments or the *Synimed* endomedullary nails.
- The surgeon and the staff at the O.R. are accountable for not respecting the warnings mentioned in these instructions for use.
- Never strike the instruments or they will lose their precision: use the sliding hammer connected to the guide to facilitate the insertion of the nail.
- Instruments must be kept away from high magnetic fields (above 1.5 Tesla).
- Decontaminate the products before removing them from the O.R.

Precaution:

These instruments have been designed so that each of their parts is used in a smooth way. If a part becomes stuck or seized, it must never be forced: unscrew it and disassemble it, then re-assemble it until it works smoothly.

Never apply excessive force while fastening the different components: doing so would produce breakage, blocking and the loss of precision of the set.

The life span of these instruments will depend on their correct use in the O.R., on their cleaning and maintenance conditions, and on their storage in a controlled humidity environment.

Do not use in the following cases:

If appropriate sterility conditions are not assured.
If you observe damages in the parts or instruments to be used in the implanting procedure.
If you observe a malfunction in any part or instrument.

Synimed shall not be liable for any damage derived from the wrong use of these instruments, or implants.

Assembling and disassembling of instruments for femoral or tibial nails:


Procedure:

Assembling:

- A** Optional: You can screw on the Guide Handle **19**, onto the body of the guide **01** if you need to perform additional maneuvers while inserting the nail.
- B** Place the femoral coupling **03** / tibial coupling **06** – depending on the bone in question- into the upper hole of the body of the guide, and insert the coupling key into the key hole.
- C** Tighten it by means of nut **04** with the ratchet wrench **20**.
- D** Place the femoral **02** / tibial **05** coupling screw into the femoral **03** /tibial **06** coupling and then screw it onto the selected nail metal core.
- E** Then screw the sliding hammer rod **14** into the threaded hole on the upper face of the body of the guide **01**. Use the combined wrench **20** on the “cube” to fasten the rod.
- F** Slide the weight **13** along the sliding hammer rod **14**.
- G** Connect handle **15** to the sliding hammer rod **14**.
- H** Slide the locking screw sleeves **23** through the holes of the guide body **01**, and fix them into the desired position by means of the adjustment knobs.
- I** Slide the soft tissue punch **21**, and the drill bit guides **22** through the locking screw sleeves **23** as needed.

Disassembling procedure:

- J** Remove the locking screws sleeves **23**, by loosening the knobs on the body of the guide **01**.
- K** Loosen the sliding hammer handle **15** by unscrewing it from the rod **14**.
- L** Remove the weight **13**.
- M** Remove the sliding hammer rod **14**, using the combined wrench **20** to unscrew it.
- N** Remove the coupling screw **02** / **05** using the combined wrench.
- O** Loosen the nut **04** by means of the combined wrench **20**.
- P** Remove the femoral **03** or tibial **06** coupling.

Decontamination and Cleaning

Synergie Ingénierie Médicale SARL has submitted the instruments items to a decontaminating and cleaning process in order to eliminate the residues produced during their manufacturing, repair or modification. They must definitely not be sterilized without having been submitted to this process, and the health care organization will be responsible for their decontamination and cleaning before sterilization.

All demountable items must carefully be disassembled before their decontamination in order to access all complex shaped areas (cavities, orifices, grooved surfaces, etc.) and ensure the action of the decontaminating product.

The operations required to comply with this requirement are described in the corresponding paragraph. When in doubt, contact Synergie Ingénierie Médicale SARL to find the appropriate information. The same precautions apply for the re-assembling procedure.

Decontamination:

Aimed at reducing the micro-organism load to limit their contaminating power.

The instruments are immersed –disassembled– into a decontaminating liquid consisting of a detergent and a disinfectant with antibacterial, antifungal and antiviral properties. Alkaline products and corrosive solutions must be avoided as they might damage the stainless steel and aluminum devices.

In order to achieve an effective decontamination, these instruments must be immersed for at least 15 minutes, but it is not advisable to exceed this time, or they might be damaged.

Once the instruments are decontaminated, they should be carefully rinsed with demineralized, de-ionized cold water. In order to increase the mechanical effectiveness of the process, it is advisable to use a non- metal brush. All these operations must be carried out before the cleaning procedure itself.

Cleaning:

This is performed through a combination of a mechanical action (rubbing, brushing), a chemical action (cleaning agent), and a thermal action (the high temperature reduces the cleaning time), which permits to eliminate dirt.

The most commonly used cleaning methods are: manual, by means of an aspersion machine, through immersion and by ultrasounds.

The choice of the cleaning method corresponds to the authorized staff in the Healthcare care facility. Do not use alkaline or corrosive products to clean

devices in stainless steel or aluminum. Special attention must be paid when cleaning instruments with complex shapes (cavities, orifices, grooved surfaces, etc.), making sure that these devices are thoroughly cleaned. Once the cleaning procedure is finished, instruments must be rinsed and dried (either manually or by means of a stove).

Certain areas of the instruments may need to be dried by means of medical compressed air.

For manual drying you must use a clean, soft, absorbing, lint free cloth.

Sterilization:

In order to prevent a re-contamination or damage to the instruments, the devices should be immediately put back in their corresponding places inside the boxes and trays. These boxes and trays must be submitted to the same decontaminating and cleaning processes as the instruments contained there-in. All the above mentioned elements must undergo a sterilization process according to the standards established by the Healthcare organization, the European Pharmacopoeia and the national regulations in force.

Synergie Ingénierie Médicale SARL recommends an autoclave sterilization at a temperature of 134 °C and a 2 bar pressure for 18 minutes. *Synimed* does not assure the preservation of the sterility of the elements which contain the instruments. According to the regulations in force, this sterility and its preservation are under the responsibility of the Healthcare organization.

Regardless of the decontamination, cleaning, and sterilization procedures implemented by the sanitary organization, it will be convenient to have them validated.

Storage, transportation and preservation:

Between two surgical procedures, the devices must be stored in a clean, dry and temperate environment/ place. They must be protected against falls, strikes and any other situation which might damage their integrity, cleanliness and/or sterility.

Their handling and transportation must be performed in such a way as to protect them from strokes or falls which might damage them. It is recommended not to store or transport boxes which are not correctly closed, and also not to pile them up as it might damage their lids and/ or seals.

If there is any doubt as to the partial or total integrity of a box, it must immediately be sent to Synergie Ingénierie Médicale SARL, in its original packaging.

Maintenance, modification and repairing:

Before returning them to the Healthcare organization,

Synergie Ingénierie Médicale SARL shall verify the operation of the instruments after having performed their corresponding maintenance, modification or repairing.

A follow up card, concerning the composition of the ancillary and the detailed technical verifications performed, shall be attached to the instruments returned to the Healthcare facility.

Once the instruments are received and before they are used in a surgical procedure, their adequate functioning must be controlled by qualified personnel authorized by the Healthcare organization.

Defective instruments should immediately be returned to Synergie Ingénierie Médicale SARL or to a firm authorized to this end. The instruments must not be repaired or modified by any person or firm not authorized by Synergie Ingénierie Médicale SARL. Only Synergie Ingénierie Médicale SARL or a firm authorized by Synergie Ingénierie Médicale SARL can perform these operations so as to guarantee the functioning and performance of the instruments in question.

Synergie Ingénierie Médicale SARL may decide not to repair an instrument, and to replace it by a new one or by one in similar conditions and working order. All verbal communication between Synergie Ingénierie Médicale SARL staff and the Healthcare organization must be subsequently confirmed in writing.

Instruments or accessories sent to Synergie Ingénierie Médicale SARL: Before being returned to *Synimed*, they must be submitted to the necessary decontaminating and cleaning processes. These operations must be recorded in the follow up card designed to this end,

which must also be sent to Synergie Ingénierie Médicale SARL.

Any exchange or transfer of these instruments between two Healthcare organizations, without the authorization of Synergie Ingénierie Médicale SARL, is strictly forbidden.

Attention:


If the Healthcare organization, the surgeon, or the personnel in the O.R. become aware that the instruments have been used in patients who are carriers of severe diseases, with high spread risks, ruled by regulations such as in the case of NCTA (non conventional transmissible agents), the said instruments will be submitted to the treatment required by law/regulations, and *Synimed* must be informed of the procedure performed.

Synimed shall not be held responsible for the actions of the Healthcare facility regarding the application of the abovementioned procedures.

Assistance and associated services:

Contact Synergie Ingénierie Médicale SARL to obtain all complementary information, technical explanations, advice or assistance regarding the performance or use of whichever instrument or accessory.

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Synergie Ingénierie Médicale S.A.R.L.