Synicem Cranioplastie

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DESCRIPTION:

SYNICEM CRANIOPLASTIE cements is self-polymerising acrylic compound with spontaneous solidification, consisting of a powder fraction formed by a methyl methacrylate polymer with a polymerisation starter, Benzoyl Peroxide, and a liquid fraction formed by Methyl Methacrylate stabilized with Hydroquinone and activated with a N,N-dimethyl p-toluidine.

When the monomer and the polymer come into contact, and are properly mixed, the N,N Dimethyl p-toluidine activates the Benzoyl Peroxide which in turn starts the polymerisation of the monomer. The polymerisation process is gradual and the mixture, fluid at first, becomes viscous and finally reaches a solid state. The reaction is exothermic, with variable maximum temperatures reaching 80°C, depending on volume and dissipation, according to thickness.

DIRECTIONS FOR USE:

SYNICEM CRANIOPLASTIE is specially used in Neurosurgery, Maxillo Facial Surgery, Othorhinolaringology and Plastic Surgery for the repair of facial and cranial defects, of traumatic, surgical or tumoral origin.

PRESENTATION:

SYNICEM CRANIOPLASTIE is supplied in a double dose pack inside a blister sterilized by ethylene oxide.

Each dose consists of:

- 1 ampoule with 17 ml of liquid sterilized by ultrafiltration in a blister sterilized by ethylene oxide.
- 1 envelope with 30 g of acrylic powder sterilized by ethylene oxide.
- 1 Sterile plastic sleeve.

CAUTIONS FOR USE:

Ensure that the temperature in the operating room does not exceed 23°C; the following periods of time (with a $\pm$ 30 sec. possible variation) apply to a 23°C cement and mixing elements room temperature. The higher the temperature, the shorter the time for working the cement and cement curing and the lower the temperature, the longer the time (see table).

The introduction of air in the cement should be avoided since bubbles may be formed in the mass. This should be avoided in the mixture (mix smoothly without lifting the spatula from the bottom of the container) as well as in the decanting or modelling.

PREPARATION AND APPLICATION:

1. Empty the blister on a sterile surface
2. Open the liquid ampoule and empty the contents into a sterile stainless steel, polythene or glass container with a capacity of at least 130 ml. The container and the spatula must be dry and at room temperature. High temperature on these elements causes premature hardening of SYNICEM CRANIOPLASTIE.
3. Open the powder envelope with sterilized scissors and empty into the liquid. The powder/liquid

Ratio is balanced in order to reach maximum material performance.

4. Mix both components until the mass reaches dough-like consistency, approximately during 30 seconds. In 2 minutes more, it can be held with the hands to mix and model it. If it sticks to the gloves, mix it a few more seconds until stickiness disappears. The material is ready to be used when it acquires certain consistency and does not stick to the gloves.
5. The plastic dough must be compressed between two flat glass surfaces (previously moisturized with sterile physiologic solution) pressing them until the desired sheet thickness is obtained. The cement mass may also be placed inside one of the sterile sleeve provided to obtain the desired thickness by rolling it with a stainless steel roller.
6. Once the sheet is obtained it must be trimmed in the shape and size required, leaving an excess edge, which is to be trimmed once more (with scissors or surgical knife) in the shape of the bone defect to be repaired.

Caution: Do not let polymerise the piece obtained “in situ” since the heat released by the acrylic when curing may injure tissues. The temperature reached is directly related to thickness (the thicker, the higher the temperature) and it can exceed 80°C.

Once polymerised, the prosthesis may be pierced with a drill, bit or punch to make the holes needed for suture or granulation. It may be positioned by means of non-reabsorbable suture material, cables, wires or titanium plate and screw systems.

If small volume or thickness defects are to be repaired, the cement can be let polymerise “in situ”, taking the care of touching to control the temperature and irrigating the material with cold physiologic solution.

The setting time at 23°C is approximately 13 - 15 minutes from the start of the mix. The setting time is also dependent of the thickness and of the quantity of resin used.

Temperature	Mixing Time	Manipulation begin	Net working time	Setting time from the start of the mix
23° C	30 sec	máx 5 min	(5 to 7) min	(13 - 15) min

PRECAUTIONS:

Monomers are volatile liquids that should be handled with special care. Prevent unnecessary contact with skin and mucosae. Carefully open ampoule avoiding splashing the contents. Pour content immediately into the mixing bowl mentioned before and add the polymer. Mix it in such a position exposure to vapours is avoided so as to prevent possible irritations.

Use during pregnancy: No data from teratogenic studies in humans are available. Product use during pregnancy or on women who may be pregnant requires a risk /benefit evaluation.

ADVERSE REACTIONS:

No allergic reactions have been reported to components of the piece set.

PHARMACEUTICAL PRECAUTIONS:

- Do not use this product after expiry date shown on the product label.
- Single use.
- DON'T REUSE. Single use product, to use for only one patient: their it reuses it can involve risks mechanic, physicochemical and/or of biological contamination.

- Handle with care and store in dark and dry place below 25° C.
- Discard if any of the sterilised containers is damaged or open.
- Do not resterilize any component.
- Do not alter powder/liquid ratio.

Formulas:

Each powder envelope contains	
Methyl Polymethacrylate	29,49 g
Benzoyl Peroxide	0,51 g

Each liquid ampoule contains	
Methyl Methacrylate	16,8 ml
N,N dimethyl p-toluidine	0,2 ml
Hydroquinone	18-20 ppm

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Synergie Ingénierie Médicale S.A.R.L.