# Bioresorbable pins

**CE 2797**

**en** This instruction for use can be downloaded as a pdf file from the website www.eifu-online.com or by flashing the QR code above, and using the UDI-DI number on the product label or the reference of this IFU. The website is accessible with the following updated browsers: Chrome, Edge, Firefox, Opera, Safari. A printed version of this instruction for use can be sent on request, without charge, under 7 calendar days by asking to eifu@teknimed.com.

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#### GENERAL INFORMATION

**Before using products, the healthcare professional must thoroughly familiarise with the safety rules mentioned in the current instructions for use, as well as the information relating to each product (description, leaflets, etc.). The corresponding information is available from qualified distributors. The healthcare professional should also be informed of the residual risk(s) associated with the product considered to be used.**

Summary of Safety and Clinical Performance is available on EUDAMED database through the link https://ec.europa.eu/tools/euamed/ searched by searching with the Basic UDI-DI mentioned above; or through https://myimplantcard.com/sscp/ until the functionality is available on EUDAMED.

#### INTENDED USERS

Medical devices should only be used in an adapted environment and by qualified healthcare professionals (surgeons) having a sound knowledge and full mastery in preparation techniques specific to such products. Preparation techniques may be acquired from qualified distributors. The healthcare professional is responsible for any complications or harmful consequences which might result from an erroneous indication or operative technique, an improper use of the equipment and/or a failure to comply with the safety rules provided in the instructions for use. Neither manufacturer nor qualified distributor can be held responsible for these complications.

#### INTENDED USE

The **Bioresorbable pins** are intended for the stabilization of metatarsal and phalangeal realignment in osteotomy procedures.

#### LIFETIME

**Bioresorbable pins** are an implantable device, intended to be totally absorbed after an average of 4 years. This period may vary depending on the physiological condition of the patient. The device is not intended to be removed and does not need any maintenance, unless a medical complication would require a surgery.

#### INDICATIONS

The **Bioresorbable pins** are indicated for the treatment of hallux valgus.

#### TARGET POPULATION

Adult patients suffering from hallux valgus.

#### PERFORMANCES AND EXPECTED CLINICAL BENEFITS

The device allows treatment of hallux valgus resulting in pain alleviation and improvement of quality of life of the patient.

#### CONTRAINDICATIONS

- Procedures other than those stated in the INDICATIONS section.
- Insufficient bone quantity or quality which can compromise correct anchoring of the pin.
- Patients susceptible to allergic reactions to the components of the device and its products of metabolism.
- Due to non-sufficient available clinical evidence, the device is contraindicated in paediatric patients, and breast-feeding or pregnant women.

#### SIDE EFFECTS AND ADVERSE EVENTS

- Complications usually encountered with bioabsorbable implants: inflammatory reaction.
- Complications usually encountered with pin: failure of the fixation, migration and protruding of the pin, pin breakage.
- Other reported side effects include: pain, swelling, wound healing issues.

Any serious incident occurring in relation to the device must be reported without delay to the manufacturer and the competent local authority where the user and/or the patient is established.

#### COMPOSITION

100% Poly(70/30; L/DL)lactide		
<b>Designation</b>	<b>Diameter</b>	<b>Reference</b>
<b>Bioresorbable pins</b>	2mm	T7100220
<b>Bioresorbable pins</b>	2.4mm	T7100224

#### STERILISATION

The product is delivered sterile under double packaging, ready to use in the operating room. The product is sterilized using Ethylene Oxide. Sterility is guaranteed until the expiry date if the packaging has not been opened or damaged.

#### RECOMMENDED ACCESSORIES

It is recommended to use **Bioresorbable pins** with the dedicated **Stainless Steel Pins** (not included).

<b>Designation</b>	<b>Diameter</b>	<b>Reference</b>
<b>Stainless Steel Pins</b>	2mm	T067120
<b>Stainless Steel Pins</b>	2.4mm	T067124

#### INSTRUCTIONS FOR USE

For the proper use of the **Bioresorbable pins**, it is essential to establish beforehand the osteotomy procedures to follow before, during, and after the operation, and also to establish the type and the location of the implant to be used.

The pin is introduced into a tunnel previously drilled out in the bone using a **Stainless Steel Pins** of the same diameter. The pin may then be cut to the appropriate length using a pair of cutting pliers.

#### WARNINGS AND PRECAUTIONS FOR USE

Read instructions for use carefully prior to use and follow preparation and handling instructions of **Bioresorbable pins**.

Ignoring the instructions for use may lead to potential undesirable effects.

For an optimum use of the device, it is essential to perform a thorough pre-operative review of the patient before the intervention in order to confirm the indication and plan the surgical technique.

Conservation

**- It is strictly forbidden to re-sterilise the product.** This product is

sterilized sterile only package has been opened or damaged.

- This device is packaged and sterilised for single use only. Do not reuse, reprocess or resterilise. Reuse, reprocessing, or resterilisation may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury or illness. Also, reprocessing or resterilisation of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another.

#### CONTRE-INDICATIONS

- Before use, check the protective packaging carefully to ensure that it has not been opened nor damaged in a way that could affect its sterility.

- Do not use this product after the expiration date printed on the package.

Use

- Examine the product visually to identify any defects such as crack or deformation. Do not implant products with defects.

- When removing the product from its packaging, be sure to follow asepsis rules.

- If packaging is unintentionally opened before use or damaged, do not use the product.

#### INFORMATION TO THE PATIENT

The patient should be informed by the healthcare professional of potential consequences of factors mentioned in the paragraphs « CONTRA-INDICATIONS » and « SIDE EFFECTS AND ADVERSE EVENTS », meaning those liable to hinder the success of the operation, as well as the possible complications which may arise. The patient should also be informed of measures to be taken to diminish possible consequences of these factors.

The patient must be informed by the healthcare professional that after implantation, this medical device does not need any precautions or specific measures regarding to environmental conditions (such as security portals in airports, induction cooking plates, etc.), or any medical examinations (such as CT-scan or MRI imaging). The medical device is indeed composed of materials that are electrically nonconductive, nonmetallic, and nonmagnetic. The medical device is thus considered MR safe (safe to be used in Magnetic Resonance environment).

#### STORAGE

Each product must be stored unopened in its original packaging at a temperature below 25°C.

#### RECOMMENDATIONS FOR DISPOSAL

Comply with local regulations in force related to handling and disposal of the products and its packaging.

#### IMPLANT CARD

As part of traceability of its devices, the manufacturer is providing for the patient an implant card on which a label specific to the device must be applied. The implant card as well as the label are located in the device box. Patients shall keep this implant card at all times in their medical records.

Instruction for completion

- **?** : Name of the patient or patient ID. To be filled by the healthcare institution.
- **[ ]** : Date of implantation. To be filled by the healthcare institution.
- **■** : Name and address of the healthcare institution. To be filled by the healthcare institution.
- "IMPLANT CARD" sticker with device information. To be affixed by the healthcare institution.

Additional information for the patient about the implanted device is available on the following website: https://myimplantcard.com.

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