

Synicem

40C

Hip Spacer
Spacer de Hanche
Espaciador de Cadera
Spaziatore d'anca
Hüft-Spacer

Head + neck length
Longueur de la tête + col
Largo cabeza + cuello
Lunghezza della testa + collo
Schenkelhals und kopf länge
48 mm

40 mm Head Diameter

Charnley Stem Type

Diamètre de la tête 40 mm

Tige type Charnley

Diámetro de Cabeza 40 mm

Vástago tipo Charnley

Diametro della testa 40 mm

Stelo tipo Charnley

40 mm Kopfdurchmesser

Charnley-Schaft

127 mm

CAUTION ALWAYS CHOOSE TO PRINT IN "ACTUAL SIZE" IN ORDER TO MAINTAIN MAGNIFICATION (NEVER FIT/SHRINK TO PRINTABLE AREA).

MAGNIFICATION / GROSSISSEMENT / INGRANDIMENTO / VERGRÖßERUNG: 15%



Synicem

48C

Hip Spacer
Spacer de Hanche
Espaciador de Cadera
Spaziatore d'anca
Hüft-Spacer

Head + neck length
Longueur de la tête + col
Largo cabeza + cuello
Lunghezza della testa + collo
Schenkelhals und kopf länge
58 mm


48 mm Head Diameter
Charnley Stem Type

Diamètre de la tête 48 mm
Tige type Charnley

Diámetro de Cabeza 48 mm
Vástago tipo Charnley

Diametro della testa 48 mm
Stelo tipo Charnley

48 mm Kopfdurchmesser
Charnley-Schaft



CAUTION ALWAYS CHOOSE TO PRINT IN "ACTUAL SIZE" IN ORDER TO MAINTAIN MAGNIFICATION (NEVER FIT/SHRINK TO PRINTABLE AREA).
MAGNIFICATION / GROSSISSEMENT / INGRANDIMIENTO / VERGRÖßERUNG: 15%

Synicem

48M

Hip Spacer
Spacer de Hanche
Espaciador de Cadera
Spaziatore d'anca
Hüft-Spacer

Head + neck length
Longueur de la tête + col
Largo cabeza + cuello
Lunghezza della testa + collo
Schenkelhals und kopf länge
58 mm


48 mm Head Diameter
Müller Stem Type

Diamètre de la tête 48 mm
Tige type Müller

Diámetro de Cabeza 48 mm
Vástago tipo Müller

Diametro della testa 48 mm
Stelo tipo Müller

48 mm Kopfdurchmesser
Müller-Schaft



CAUTION ALWAYS CHOOSE TO PRINT IN "ACTUAL SIZE" IN ORDER TO MAINTAIN MAGNIFICATION (NEVER FIT/SHRINK TO PRINTABLE AREA).

MAGNIFICATION / GROSSISSEMENT / INGRANDIMENTO / VERGRÖßERUNG: 15%



Synicem

56C

Hip Spacer
Spacer de Hanche
Espaciador de Cadera
Spaziatore d'anca
Hüft-Spacer

Head + neck lenght
Longueur de la tête + col
Largo cabeza + cuello
Lunghezza della testa + collo
Schenkelhals und kopf länge
65 mm

56 mm Head Diameter
Charnley Stem Type

Diamètre de la tête 56 mm
Tige type Charnley

Diámetro de Cabeza 56 mm
Vástago tipo Charnley

Diametro della testa 56 mm
Stelo tipo Charnley

56 mm Kopfdurchmesser
Charnley-Schaft

130 mm

CAUTION ALWAYS CHOOSE TO PRINT IN "ACTUAL SIZE" IN ORDER TO MAINTAIN MAGNIFICATION (NEVER FIT/SHRINK TO PRINTABLE AREA).

MAGNIFICATION / GROSSISSEMENT / INGRANDIMENTO / VERGRÖßERUNG: 15%



Synicem

56M

Hip Spacer
Spacer de Hanche
Espaciador de Cadera
Spaziatore d'anca
Hüft-Spacer

Head + neck lenght
Longueur de la tête + col
Largo cabeza + cuello
Lunghezza della testa + collo
Schenkelhals und kopf länge
65 mm

56 mm Head Diameter
Müller Stem Type

Diamètre de la tête 56 mm
Tige type Müller

Diámetro de Cabeza 56 mm
Vástago tipo Müller

Diametro della testa 56 mm
Stelo tipo Müller

56 mm Kopfdurchmesser
Müller-Schaft

